

FILLING HANDLER - General Standard Work & Information

Line 10 Batch Change Standard Work	LPF 1 & 2	p. 2-4
	LPF 3 & 4	p. 5-7
	Tank Room	p. 8-9
	Tank Mixing	p. 10
Line 11 Batch Change Standard Work	Lead / Tank Room Handler	p. 11
	Night Lead/Tank Room Handler	p. 12
	LPF 5, 6 and 7 Handler	p. 13
	RPF 1 and 2 Handler	p. 14
	DAY Shift LPF/RPF Operator & Day Shift Manual Belt Operators	p. 15
	Manual Belt Handler	p. 16
	Tank Room Prep	p. 17
Line 11 Batch Change with Sanitization Standard Work	Lead / Tank Room Handler	p. 18
	Night Lead/Tank Room Handler	p. 19
	LPF 5, 6 and 7 Handler	p. 20
	RPF 1 and 2 Handler	p. 21
	Day Shift LPF/RPF Operator & DAY Shift Manual Belt Operators	p. 22
	Night Manual Belt Handler	p. 23
	Manual Belt Handler	p. 24
	Tank Room Prep	p. 25
Line 10 Tool Crib Std. Work		p. 26
Line 11 Tool Crib Std. Work		p. 27-33
Tool Re-order Sheet		p. 34
Tool Shift Verification		p. 35
Line 10 Parts Kanban		p. 36
Enhanced Closure Go/No Go Gauge Standard Work		p. 37-38
Solution Pump Lockout Std. Work		p. 39
Rework % Reduction		p. 40
Cyclohex Dispensing		p. 41
Line 11 Piping Verification Talley Sheet		p. 42-43
Line 10 6S Walk-around		p. 44
Line 11 6S Scorecard		p. 45
FOB Shift Checks		p. 46
LPF and RPF – Maintenance, Flush & Sampling Notification		p. 47

Line 10 “Fill Side” Batch Change Standard Work

LPF01

External

- Verify imprint station is setup with old batch BOM, setup sheet, imprint board and / or calibrated ruler.
- Get heat resistant gloves, drill and drill extension for imprint station.
- Get type box for LPF once the new batch BOM / setup sheet has been locked in box.
- Get barcode cords plugged in and ready (only if changing codes).
- Get manifolds and assure tubing is in good working condition.
- Keep rework bin cleaned out.
- The handler that performs the last tank change of the batch will call the Chem Lab and notify them that we are in our last tank, (this is an Opti-Lab requirement).
- The handler that performs the last quality check of the batch will close old batch out in SPC system.
- Perform FOB checks and document in fill / pack.

Internal

- Cycle stop machine, turn accumulator off, use clear machine function and let machine clear out until last pattern of bags are above nozzles or at end of barcode station. Cycle stop machine, use clear machine function in Imprint Mode. Cycle stop / E-stop LPF after removing the last pattern of unfilled printed bags from end of LPF conveyor which will be used as EOB imprints.
- Handlers obtain EOB imprints, check and sign.
- Place BOM / setup sheets & EOB imprints in old batch folder (after last imprints are pulled and signed).
- Move the fill station downhill away from barcode station if applicable, to get LPF ready to flush.
- Remove old type from the die plate and lock it up in the type removed side of the type box. Once the old type is locked up, unlock the type to be installed side of the box.
- Proceed to the unload station between LPF02 & LPF03. Line Clearance can begin when all handlers are present. Maintain proper interval during clearance. Perform line clearance on LPF02 & LPF01. Document clearance in Fill / Pack.
- As soon as the Clearance is complete, disengage E-Stops and push machine reset to clear all fault messages. On the HMI screen “select” clear machine, mode of operation-automatic, accumulator on, imprint mode enabled than “select” clear entire machine. Run machine (8 cycles) until printed, unfilled bags are one cycle past the fill station. Cycle stop the machine and depress the E-Stop. Remove the set of bags that are past the fill station (close clamps) for FOB Imprints, check and sign if acceptable.
- Install manifold / verify air inject air lines are removed (manual nozzle also).
- Perform FOB flush if not previously performed (for FOB flush, machines must be flushed in sequence). Ensure flush gets documented in fill/pack.
- Remove flush manifold and hook up air inject air lines (manual nozzle also).
- When all machines have completed the flush, start LPF and pull FOB samples while belt is running, place in sample cart.

Operators

- Operators help handler obtain EOB imprints.
- Operators remove used hot stamp and barcode foil (used hot stamp foil is placed on foil cart and used barcode foil is placed in trash).
- Verify reject bins are empty and rejects are cut.
- Take filling and pouching trash to appropriate pouching side door.
- Perform Zone Clearance.

FOR INFORMATION ONLY – MAY BE DISCARDED **DOCUMENT OWNER: N. BUCHANAN**

1Copy – LPF01 Attached to clipboard

Line 10 “Fill Side” Batch Change Standard Work

LPF02

External

- Verify imprint station is setup with old batch BOM, setup sheet, imprint board and / or calibrated ruler.
- Get heat resistant gloves, drill and drill extension for imprint station.
- Get type box for LPF once the new batch BOM / setup sheet has been locked in box.
- Get barcode cords plugged in and ready (only if changing codes).
- Get manifolds and assure tubing is in good working condition.
- Get sample cart for LPF01 / LPF02.
- Keep rework bin cleaned out.
- The handler that performs the last tank change of the batch will call the Chem Lab and notify them that we are in our last tank, (this is an Opti-Lab requirement).
- The handler that performs the last quality check of the batch will close old batch out in SPC system.
- Perform FOB checks and document in fill / pack.

Internal

- Cycle stop machine, turn accumulator off, use clear machine function and let machine clear out until last pattern is above nozzles or at end of barcode station. Cycle stop machine, use clear machine function in Imprint Mode. Cycle stop / E-stop LPF after removing the last pattern of unfilled printed bags from end of LPF conveyor which will be used as EOB imprints.
- Handlers obtain EOB imprints, check and sign.
- Place BOM / setup sheets & EOB imprints in old batch folder (after last imprints are pulled and signed).
- Move fill station downhill away from barcode station if applicable, to get machine ready to flush.
- Remove old type from the die plate and lock it up in the type removed side of the type box. Once the old type is locked up, unlock the type to be installed side of the box.
- Proceed to the unload station between LPF02 & LPF03. Line Clearance can begin when all handlers are present. Maintain proper interval during clearance. Perform line clearance on LPF02 & LPF01. Document clearance in Fill / Pack.
- As soon as the Clearance is complete, disengage E-Stops and push machine reset to clear all fault messages. On the HMI screen “select” clear machine, mode of operation-automatic, accumulator on, imprint mode enabled than “select” clear entire machine. Run machine (8 cycles) until printed, unfilled bags are one cycle past the fill station. Cycle stop the machine and depress the E-Stop. Remove the set of bags that are past the fill station (close clamps) for FOB Imprints, check and sign if acceptable.
- Install manifold / verify air inject air lines are removed.
- Perform FOB flush if not previously performed (for FOB flush, machines must be flushed in sequence). Ensure flush gets documented in fill/pack.
- Remove flush manifold and hook up air inject air lines.
- When all machines have completed the flush, start LPF and pull FOB samples while belt is running, place in sample cart.

Operators

- Operators help handler obtain EOB imprints.
- Operators remove used hot stamp and barcode foil (used hot stamp foil is placed on foil cart and used barcode foil is placed in trash).
- Verify reject bins are empty and rejects are cut.
- Take filling and pouching trash to appropriate pouching side door.
- Perform Zone Clearance.

FOR INFORMATION ONLY – MAY BE DISCARDED

DOCUMENT OWNER: N. BUCHANAN

1Copy – LPF02 Attached to clipboard

Line 10 “Fill Side” Batch Change Standard Work

LPF01 / LPF02 Helper

External

- Get Walkie-Talkie’s and drill extensions from office, have available for batch change.
- Get new cleaning papers and ticket pouches ready for new batch.
- Get paperwork for new batch from dry side.
- Get type boxes ready with new BOM and setup sheets. Verify type is secure in type retainers. Lock paperwork up inside the type to be installed side of the box.
- Lock any additional paperwork up inside the big die plate cabinet.
- Verify that the SPC has been closed out after the last quality check of the batch is completed.
- Document all “partial” closure totes, barcode foil, hot stamp foil and bag closets in the material consumption portion of fill/pack.
- Get old batch folder ready. Verify that the following items are in/or ready to be placed in the batch folder. (Previous days cleaning forms, major maintenance, bag, closure and foil tickets, FOB imprints and any adjustment imprints).
- If Opti-Lab is operational, verify that the Chem Lab was notified when the last tank of the batch was started.
- If applicable, ensure new air inject for manual nozzle is ready to be installed at batch change.
- Clean top of LPF01/02 sample cabinet with alcohol. Place the tagged bags that are going to be used for EOB samples on top of the cabinet.

Internal

- Use Walkie-Talkies to communicate with person in tank room, (i.e., samples pulled, clear machines out, put flush manifolds on, flush machines, etc).
- Enter any remaining bag scrap into fill/pack and help remove scrap / trash from wet side.
- Spray gloves and manual nozzle with alcohol. Pull EOB samples after all machines have cleared out. Verify samples in fill/pack prior to filter time expiring if applicable.
- Prepare manual nozzle for flushing, (switch nozzle to off position, remove air inject line, change air inject if applicable and place flush pipe under nozzle).
- Remove the air inject lines from nozzles on LPF02.
- Remove the air inject lines from nozzles on LPF01.
- Place EOB imprints and any other batch specific documents in the batch folder. Hand the batch folder out to the dry side.
- Start with LPF02, insert the new batch type retainers into the die plate, rethread black foil, go to HMI screen and “select” batch change selection, “select” enter new batch mode than “select” OK, (verify that the counters reset to zero).
- Reprogram the barcode information only if the Code is changing.
- Go to LPF01 and perform the same steps as stated for LPF02.
- As soon as the clearance on LPF02 is complete, push machine reset to clear all fault messages. On the HMI screen “select” clear machine, mode of operation-automatic, accumulator on, imprint mode enabled than “select” clear entire machine. Run machine (8 cycles) until printed, unfilled bags are one cycle past the fill station. Cycle stop the machine and depress the E-stop. Remove the set of bags (a pattern of 5) that are past the fill station (close clamps) for FOB imprints, check and sign.
- Help install / remove flush manifolds and hookup air inject lines as needed.
- Go to LPF01 and assist with the imprint process as listed above for LPF02.
- Be available to perform 2nd verification on FOB imprints as needed.

Line 10 “Fill Side” Batch Change Standard Work

LPF03

External

- Verify imprint station is setup with old batch BOM, setup sheet, imprint board and / or calibrated ruler.
- Get heat resistant gloves, drill and drill extension for imprint station.
- Get type box for LPF once the new batch BOM / setup sheet has been locked in box.
- Get barcode cords plugged in and ready (only if changing codes).
- Get manifolds and assure tubing is in good working condition.
- Get sample cart for LPF03 / LPF04.
- Keep rework bin cleaned out.
- The handler that performs the last tank change of the batch will call the Chem Lab and notify them that we are in our last tank, (this is an Opti-Lab requirement).
- The handler that performs the last quality check of the batch will close old batch out in SPC system.
- Perform FOB checks and document in fill / pack.

Internal

- Cycle stop machine, turn accumulator off, use clear machine function and let machine clear out until last pattern of bags are above nozzles or at end of barcode station. Cycle stop machine, use clear machine function in imprint mode. Cycle stop / E-stop LPF after removing the last pattern of unfilled printed bags from end of LPF conveyor which will be used as EOB imprints.
- Handlers obtain EOB imprints, check and sign.
- Place BOM / setup sheets & EOB imprints in old batch folder (after last imprints are pulled and signed).
- Move the fill station downhill away from barcode station if applicable, to get LPF ready to flush.
- Remove old type from the die plate and lock it up in the type removed side of the type box. Once the old type is locked up, unlock the type to be installed side of the box.
- Proceed to the unload station between LPF02 & LPF03. Line Clearance can begin when all handlers are present. Maintain proper interval during clearance. Perform line clearance on LPF03 & LPF04. Document clearance in Fill / Pack.
- As soon as the Clearance is complete, disengage E-Stops and push machine reset to clear all fault messages. On the HMI screen “select” clear machine, mode of operation-automatic, accumulator on, imprint mode enabled than “select” clear entire machine. Run machine (8 cycles) until printed, unfilled bags are one cycle past the fill station. Cycle stop the machine and depress the E-stop. Remove the set of bags that are past the fill station (close clamps) for FOB Imprints, check and sign if acceptable.
- Install manifold / verify air inject air lines are removed.
- Perform FOB flush if not previously performed (for FOB flush, machines must be flushed in sequence). Ensure flush gets documented in fill/pack.
- Remove flush manifold and hook up air inject air lines.
- When all machines have completed the flush, start LPF and pull FOB samples while belt is running, place in sample cart.

Operators

- Operators help handler obtain EOB imprints.
- Operators remove used hot stamp and barcode foil (used hot stamp foil is placed on foil cart and used barcode foil is placed in trash).
- Verify reject bins are empty and rejects are cut.
- Take filling and pouching trash to appropriate pouching side door.
- Perform Zone Clearance.

FOR INFORMATION ONLY – MAY BE DISCARDED

DOCUMENT OWNER: N. BUCHANAN

1Copy – LPF03 Attached to clipboard

Line 10 “Fill Side” Batch Change Standard Work

LPF04

External

- Verify imprint station is setup with old batch BOM, setup sheet, imprint board and / or calibrated ruler.
- Get heat resistant gloves, drill and drill extension for imprint station.
- Get type box for LPF once the new batch BOM / setup sheet has been locked in box.
- Get barcode cords plugged in and ready (only if changing codes).
- Get manifolds and assure tubing is in good working condition.
- Keep rework bin cleaned out.
- The Handler that performs the last tank change of the batch will call the Chem Lab and notify them that we are in our last tank, (this is an Opti-Lab requirement).
- The Handler that performs the last quality check of the batch will close old batch out in SPC system.
- Perform FOB checks and document in fill / pack.

Internal

- Cycle stop machine, turn accumulator off, use clear machine function and let machine clear out until last pattern of bags are above nozzles or at end of barcode station. Cycle stop machine, use clear machine function in imprint mode. Cycle stop / E-stop LPF after removing the last pattern of unfilled printed bags from end of LPF conveyor which will be used as EOB imprints.
- Handlers obtain EOB imprints, check and sign.
- Place BOM / setup sheets & EOB imprints in old batch folder (after last imprints are pulled and signed).
- Move the fill station downhill away from barcode station if applicable, to get LPF ready to flush.
- Remove old type from the die plate and lock it up in the type removed side of the type box. Once the old type is locked up, unlock the type to be installed side of the box.
- Proceed to the unload station between LPF02 & LPF03. Line Clearance can begin when all handlers are present. Maintain proper interval during clearance. Perform line clearance on LPF03 & LPF04. Document clearance in Fill / Pack.
- As soon as the Clearance is complete, disengage E-Stops and push machine reset to clear all fault messages. On the HMI screen “select” clear machine, mode of operation-automatic, accumulator on, imprint mode enabled than “select” clear entire machine. Run machine (8 cycles) until printed, unfilled bags are one cycle past the fill station. Cycle stop the machine and depress the E-stop. Remove the set of bags that are past the fill station (close clamps) for FOB imprints, check and sign if acceptable.
- Install manifold / verify air inject air lines are removed.
- Perform FOB flush if not previously performed (for FOB flush, machines must be flushed in sequence). Ensure flush gets documented in fill/pack.
- Remove flush manifold and hook up air inject air lines.
- When all machines have completed the flush, start LPF and pull FOB samples while belt is running, place in sample cart.

Operators

- Operators help handler obtain EOB imprints.
- Operators remove used hot stamp and barcode foil (used hot stamp foil is placed on foil cart and used barcode foil is placed in trash).
- Verify reject bins are empty and rejects are cut.
- Take filling and pouching trash to appropriate pouching side door.
- Perform Zone Clearance.

FOR INFORMATION ONLY – MAY BE DISCARDED **DOCUMENT OWNER: N. BUCHANAN**

1Copy – LPF03 Attached to clipboard

Line 10 “Fill Side” Batch Change Standard Work

LPF03 / LPF04 Helper

Internal

- Turn batch change light on (located in material exchange room) when we start filling out of the surge tank.
- Take old batch boards down / put new batch boards up. Place old batch boards in batch folder.
- Enter any remaining bag scrap into fill/pack and help remove scrap / trash from wet side.
- Change over fill/pack between LPF03 & LPF04 to new batch.
- Verify Group Clearances are being performed.
- Enter verifying Group Clearances into fill/pack under new batch.
- Remove the air inject lines from nozzles on LPF04.
- Push used foil cart to dirty side of material exchange room floor.
- Remove the air inject lines from nozzles on LPF03.
- Start with LPF03, insert the new batch type retainers into the die plate, rethread black foil, go to HMI screen and “select” batch change selection, “select” enter new batch mode than “select” OK, (verify that the counters reset to zero).
- Reprogram the barcode information only if the Code is changing.
- Go to LPF04 and perform the same steps as stated for LPF02.
- As soon as the clearance on LPF03 is complete, push machine reset to clear all fault messages. On the HMI screen “select” clear machine, mode of operation-automatic, accumulator on, imprint mode enabled than “select” clear entire machine. Run machine (8 cycles) until printed, unfilled bags are one cycle past the fill station. Cycle stop the machine and depress the E-stop. Remove the set of bags (a pattern of 5) that are past the fill station (close clamps) for FOB imprints, check and sign.
- Help install / remove flush manifolds and hookup air inject lines as needed.
- Go to LPF04 and assist with the imprint process as listed above for LPF02.
- Be available to perform 2nd verification on FOB imprints as needed.
- Turn batch change light off when first bags from new batch reach the pouching side.

FOR INFORMATION ONLY – MAY BE DISCARDED DOCUMENT OWNER: N. BUCHANAN
1Copy – LPF03 Attached to clipboard

Line 10 “Fill Side” Batch Change Standard Work

Tank Room

External

- Verify the Walkie-Talkie’s are working and can be used to communicate between the Tank Room and LPF Room. If Walkie-Talkie’s are not available, Face-to-Face Communications must be conducted.

Prior to Batch Change

(Filter Change)

- Put pre-filter on water and verify all filter tags match paperwork (per SOP).
- Place check mark on old filter tags to indicate that they are USED FILTERS.
- Get measuring Rod from Vat, (# on Rod should match # on surge tank) rinse it off and place it on the Filter Cart.
- Approximately twenty minutes prior to batch change, take the new filters off water and stage them where they will be installed.

(Filter Test)

- Set up filter test site as instructed in the filter certification training.
- Get measuring rod from Vat, (# on Rod should match # on surge tank) rinse it off and place it on the filter cart.

After Batch Change

- Go inside and enter the appropriate flushes, tank hook up and filters into fill pack.
- Fill out Pre-production Solution Transmission System Cleaning Record (NC-FL-029) and give to dry side along with the Truck cards.
- Update the filter expiration time and date per SOP, located inside LPF room above the desk.
- Place used Filters on water.
- Rinse off Measuring Rod and place it back in Vat.
- Clean up Filter Test Site.

Internal

- Turn off the Supply Pump and Close Tank Valve (if dumping solution).
- Close and Lock Header Box than uncap pipe leading into Header Box.
- Notify inside Handler when you start running the solution out of the Surge Tank.
- Go to mixing and sign for unhooking the tank. If dumping solution from mix tank, use Radar or measuring Rod to determine the amount of solution remaining in the mix tank.
- Pick up new batch folder and return to the tank room. (If batch folder is not ready at this time, pick up prior to printing addressograph ticket).

Internal Cont.

- If dumping solution, open tank and drain. If dumping 5,000 Liters or more, monitor to make sure solution does not back up into fill rooms.
- Monitor solution level in surge tank. When the solution reaches the holes in the side of the bell, notify the inside to clear out the LPF's. If walkie-talkie is not available, notify the inside to clear out when solution reaches the top of the bell.
- If testing filters, take off cartridge filter and put on water per SOP. If replacing filters, remove the old prefilter and cartridge filter and install new filters as long as you are not hooked to the surge tank.
- When samples have been pulled, turn off discharge pumps and verify the amount of solution left in the surge tank and drain.
- Go inside and enter tank unhook and EOB solution loss in fill pack (mix tank, surge tank and piping).
- Enter any remaining scrap if necessary than change fill/pack over to new batch.
- Once on new batch, select fill point setup than select new tank which will print your addressograph ticket.
- Return to the tank room, Verify surge tank is empty, hang measuring rod inside surge tank and close surge tank valve.
- Test and reinstall cartridge filter or remove old multiplate filters and install new multiplate filters. Reconnect cartridge filter to surge tank (make sure all hoke valves are open).
- Have second person verify tank hook up. Perform tank valve flush and cap off tank valve per SOP.
- Open tank valve and header box. Record time in use on prefilter and cartridge filter tags per SOP.
- Turn supply pump on and do surge tank flush per SOP. Drain surge tank and remove measuring rod.
- Verify surge tank is empty, close and lock lid, close surge tank valve and start refilling surge tank.
- Record time in use on new multiplate filter tags per SOP. Once solution is coming through the hoke valves, close valves half way before turning on discharge pumps.
- Once solution starts coming through hoke valves notify inside to put manifolds on LPF's.
- Once surge tank is about half full notify inside to flush machines.

FOR INFORMATION ONLY – MAY BE DISCARDED DOCUMENT OWNER: N. BUCHANAN

1Copy – LPF01 Attached to clipboard

Batch Change Tank Mixing Standard Work for L10 Night Shift

1. Determine how fast the line is running by taking the starting volume of the current tank, allowing the tank run for 10 minutes and then recording end volume of that tank.
2. Subtract the starting volume from the end volume then divide 10. This will give you the speed at which the line is running. (Example: $(20,000-17,000) \text{ L} / 10 \text{ min} = 300 \text{ L/min}$)
3. Repeat this process several times to determine the average speed the line is running.
4. Average times are typically about 300-350 L/min.
5. Divide the tank volumes by this number to determine how many minutes the current tanks mixed will last. (Example: $39000/300 = 130$ minutes).
6. Determine what next tank volumes need to be mixed in order to allow the second to last tank of the batch to be completed by 0030.
7. The second to last tank being empty at 0030 will give mixing adequate time to finish filling and mixing the first tank of the next batch (It takes mixing at least 1.5 hours to mix the first tank). Batch Change normally begins between 02:00 and 02:30.
8. Recalculate the next tank volumes while second tank is being run to determine if new volumes are needed due to downtime. The third/fourth tank may need to be adjusted if the line has been down long enough.
9. Call mixing when the last tank change is performed or when the tank is completely empty after dumping (This is referring to the second to last tank of the current batch).
10. Remember that the Cafeteria Service lasts from 01:00-03:00; Batch Change has to begin early enough before 03:00 to give the operators a chance to eat.
11. Excessive downtime can affect when the second to last tank is empty. Be sure to have the tank empty and ready for mixing by 00:30; 01:00 at the latest. If more than 4200 L are to be dumped, call wastewater (776-6262) to determine if that amount of solution waste treatment can handle.

FOR INFORMATION ONLY – MAY BE DISCARDED

DOCUMENT OWNER: N. BUCHANAN

1Copy – LPF01 Attached to clipboard

LEAD/TANKROOM HANDLER

Batch Change Standard Work

Time	Task to be Performed
Prior to EOB	Pick up set up sheet and batch folder information.
Prior to EOB	Sign in new dies. Verify Type retainer stamps, sort paperwork, batch documentation. Place barcode setup sheets and BOM with each type box.
Prior to EOB	Clear manual printer 2199.
Prior to EOB	Sign imprint for 2199 (Printer setup prior to EOB).
Prior to EOB	Perform tank unhook to begin running out of Surge tank to allow surge to be mostly consumed by machine clear out time (18:53).
Prior to EOB	Document unhook in Fill/pack system.
Prior to EOB	Monitor Surge tank until time for machines to be cleared out/ Do not let surge tank run dry.
At EOB	Signal line to clear out (As close to 18:53 as possible).
At EOB	Once production has stopped, pull EOB samples.
At EOB	Once samples are obtained and documented, measure surge tank heel and open to drain.
At EOB	Document Surge tank heel.
At EOB	Measure mix tank heel/return key to mixing and open to drain.
At EOB	Document Mix tank heel.

Line 11 - Night Lead/Tank Room Handler

Batch Change Standard Work

(Fill Side-All Night Shift Supervisors/Asst. Supervisors/PT's/ Handlers must be garbed and on the line at 18:55)

Time	Task to be Performed
Prior to EOB	Verify release of new batch tank, select tank, perform fill point setup and print sample tags for new batch.
At FOB	Tag sample bags to be used for SP1 (tank room samples).
At FOB	Tank Room Handler and additional handler must verify the Tank Hook Up.
At FOB	Once old batch unhook has been performed and old filters removed: <ul style="list-style-type: none"> • Perform tank valve flush • Install cap at header box • Flush SP1 and pull SP1 samples
At FOB	Once new filters are installed prior to the surge and surge tank has been drained from old solution then: <ul style="list-style-type: none"> • Install surge tank measuring rod • Perform flush through Filter #1 / Filter #2 and document • Perform appropriate surge tank flush per Filling 3H and drain
At FOB	Verify Filters #3 and #4 are installed downstream of surge tank.
At FOB	Begin filling surge tank with solution for FOB flush.
At FOB	Notify Machine Handlers to install flush manifolds.
At FOB	Once surge is filled flush through Filter #3 / Filter #4 and document.
At FOB	Once surge tank has filled up then signal line to begin FOB flush per Fill 4 Addendum.
At FOB	Assure entire line has flushed the required amount of solution per Filling 4 Addendum.

Line 11 - LPF 5, 6 and 7 Handler Batch Change Standard Work

(Fill Side-All Night Shift Supervisors/Asst. Supervisors/PT's/ Handlers must be garbed and on the line at 18:55)

	Time	Task to be Performed
DAY	Prior to EOB	Have manifolds at machines, inspect for missing/ damaged tubing and correct (prior to 1853).
DAY	Prior to EOB	Get sample cart. Stage at machine (prior to 1853).
DAY	Prior to EOB	Get new batch type box and stage at the machine (prior to 1853) along with die gloves and charged drill with extension/chuck. Attain a die box key.
DAY	Prior to EOB	Get BOM, Barcode setup sheet, Barcode scanner and 'First Imprints' from the batch that is ending for imprint verification.
DAY	At EOB	To clear out machine: <ul style="list-style-type: none"> • Disable accumulator. • Clear machine button. • Once last pattern passes fill station then re-enable the accumulator to allow another bag pattern to load. • Once this pattern reaches the fill station (all other units in the machine should have been offloaded), then cycle stop and e-stop. • Assure only 1 complete pattern is pulled at a time.
DAY	At EOB	Sign last imprints.
DAY	At EOB	Place signed imprints (Last and First), BOM and Barcode setup sheet in red batch folder.
DAY	At EOB	Remove old type and install new type.
Night	FOB	Begin line clearance as soon as old type is removed and locked in box. (Should be no later than 19:00).
DAY	At EOB	Change over Fill/Pack to new batch.
DAY	At EOB	Verify and Document operator zone clearance.
DAY	At EOB	Verify previous batch SPC is closed out.
Night	FOB	Pull First Imprints. Sign Imprint and Imprint verification form.
Night	FOB	After acceptable first imprint is obtained or when Tank Room Handler announces- install manifold on machine and crack open drain valve on manifold, verify air inject lines removed.
Night	FOB	Tank Room Handler will signal to FOB flush per Filling 4 Addendum.
Night	FOB	After flush complete close hoke valve, remove manifold and reconnect air injector lines.
Night	FOB	If imprints not completed at this time continue with First Imprint process.
Night	FOB	If Imprint process complete and entire line has completed FOB flushes then perform FOB sampling.

Line 11 - RPF 1 and 2 Handler Batch Change Standard Work

(Fill Side-All Night Shift Supervisors/Asst. Supervisors/PT's/ Handlers must be garbed and on the line at 18:55)

	Time	Task to be Performed
DAY	Prior to EOB	Have manifolds at machines, inspect for missing/ damaged parts and correct (prior to 1853).
DAY	Prior to EOB	Get sample cart. Stage at machine (prior to 1853).
DAY	Prior to EOB	Get new batch type box and stage at the machine (prior to 1853) along with die gloves and charged drill with extension/chuck. Attain a die box key.
DAY	Prior to EOB	Get BOM, Barcode setup sheet, Barcode scanner and 'First Imprints' from the batch that is ending for imprint verification.
DAY	At EOB	To clear out machine: <ul style="list-style-type: none"> Assure only 1 complete pattern is pulled at a time.
DAY	At EOB	Sign last imprints.
DAY	At EOB	Place signed imprints (Last and First), BOM and Barcode setup sheet in red batch folder.
DAY	At EOB	Remove old type and install new type.
Night	FOB	Begin line clearance as soon as old type is removed and locked in box. (Should be no later than 19:00).
DAY	At EOB	Change over Fill/Pack to new batch.
DAY	At EOB	Verify and Document operator zone clearance.
DAY	At EOB	Verify previous batch SPC is closed out.
Night	FOB	Pull First Imprints. Sign Imprint and Imprint verification form.
Night	FOB	After acceptable first imprint is obtained or when Tank Room Handler announces- install manifolds on machine and crack open drain valve on fill manifold, verify air inject lines removed.
Night	FOB	Tank Room Handler will signal to FOB flush per Filling 4 Addendum.
Night	FOB	After flush complete close hoke valve, remove manifolds and reconnect air injector lines.
Night	FOB	If imprints not completed at this time continue with First Imprint process.
Night	FOB	If Imprint process complete and entire line has completed FOB flushes then perform FOB sampling.

Line 11 - DAY Shift LPF/RPF Operator

Batch Change Standard Work

Time	Task to be Performed
At EOB	Perform zone clearance of assigned machine/area. Remove any units and place in discard bin and remove batch board sheets and discard into trash or place in batch folder.
At EOB	Cut all units in the discard bin. Communicate the number of units cut to the material handler for documentation.
At EOB	Remove cut scrap/trash and place out the door to pouch side at A or E belt(Once door bar in place, door cannot be opened to pass out scrap).
At EOB	Verify handler has completed imprints then remove both the barcode foil and the black foil from the machine/printer.
At EOB	Remove air inject lines.

Line 11 - DAY Shift Manual Belt Operators

Batch Change Standard Work

Prior to EOB	Once printed bag count is achieved clear out printer. Leave last pattern of printed bags for handler's imprints. Advance and remove the black and white used foil. Perform clearance of printer area removed all printed bags and used foil from the area.
At EOB	Remove cyclohex wick applicators from nozzle area and place on rack. (Leave 1 applicator for sample nozzle).
At EOB	Disconnect air injector lines.
At EOB	Perform clearance of nozzle area and remove all trash and bag scrap through E belt door.

Line 11 - Manual Belt Handler Batch Change Standard Work

(Fill Side-All Night Shift Supervisors/Asst. Supervisors/PT's/ Handlers must be garbed and on the line at 18:55)

	Time	Task to be Performed
DAY	Prior to EOB	Determine bag lead needed.
DAY	Prior to EOB	Get sample cart. Stage at machine (prior to 1800).
DAY	Prior to EOB	Get new batch type box and stage at the machine (prior to 1800) along with die gloves and charged drill with extension/chuck. Attain a die box key.
DAY	Prior to EOB	Get BOM, Barcode setup sheet, Barcode scanner and 'First Imprints' from the batch that is ending for imprint verification (prior to 1800).
DAY	Prior to EOB	Sign last imprints.
DAY	Prior to EOB	Place signed imprints (Last and First), BOM and Barcode setup sheet in red batch folder.
DAY	Prior to EOB	Turn on Batch Change Light.
DAY	Prior to EOB	Chain off printer area.
DAY	Prior to EOB	Remove old type and install new type.
DAY	Prior to EOB	Clear 2199 printer area.
DAY	Prior to EOB	Document clearance and Print area operator zone clearance.
DAY	Prior to EOB	Pull first imprints from 2199. Verify/ sign imprint and imprint verification sheet.
DAY	At EOB	Once all production has ended pull EOB samples.
DAY	At EOB	Verify imprints present and all documents in batch folder. Place batch fold on F belt for pouching side to pick up.
DAY	At EOB	Change over Fill/Pack to new batch.
DAY	At EOB	Verify previous batch SPC is closed out.
Night	Prior to EOB	Get all material partial and cyclohex information and document materials into new batch fill/pack.
Night	At FOB	Perform FOB flush per Filling 4 Addendum when signaled by Tank Room Handler.
Night	At FOB	Printing on Multi-up may begin.
Night	At FOB	Pull FOB samples from manual nozzles.

Line 11 - Tank Room Prep Batch Change Standard Work

(Prep-All Night Shift Supervisors/Asst. Supervisors/PT's/ Handlers must be garbed and on the line at 18:50. Begins when last tank is unhooked and running out of surge)

	Time	Task to be Performed
DAY	Prior to Tank Unhook	Have new batch filters DW flushed with whirl packs on ends and ready staged on cart or by their location.
DAY	After Tank Unhook	Remove Filter #1 and Filter #2. Place whirl pack on the inlet piping to the surge tank. Test removed Filters.
Night	After Tank Unhook	As old filters are removed, install new Filter #1 and #2. <u>Do not make the final connection</u> to the surge tank until production has stopped from the previous batch and the surge tank heel has been measured and drained.
Night	At EOB	Day shift Tank Room Handler will shut off pumps, measure and open to drain. Remaining filters may be removed from the STS at this time.
Night	At FOB	Complete the final connection between Filter #2 and surge tank.
Night	At FOB	Install Filter #3 and Filter #4 on outlet side (inlet still covered with whirl packs), then complete the connection for the inlet side once the surge tank flush has been completed and drained by the Night shift Tank Room Handler.
Night	At FOB	Test Filters from the previous batch.

LEAD/TANKROOM HANDLER
Batch Change with Sanitization Standard Work

Time	Task to be Performed
Prior to EOB	Connect all DW hoses and attach to the DW water drops for sanitization.
Prior to EOB	Pick up set up sheet and batch folder information.
Prior to EOB	Sign in new dies. Verify Type retainers, stamps, sort paperwork, batch documentation. Place barcode setup sheets and BOM with each type box.
Prior to EOB	Clear manual printer 2199.
Prior to EOB	Sign imprint for 2199 (Printer setup prior to EOB).
Prior to EOB	Perform tank unhook to begin running out of Surge tank to allow surge to be mostly consumed by machine clear out time (18:00).
Prior to EOB	Document unhook in Fill/pack system.
At EOB	Monitor Surge tank until time for machines to be cleared out/ Do not let surge tank run dry.
At EOB	Signal line to clear out (As close to 18:00 as possible).
At EOB	Once production has stopped, pull EOB samples.
At EOB	Once samples are obtained and documented, measure surge tank heel and open to drain.
At EOB	Document Surge tank heel.
At EOB	Measure mix tank heel/obtain both keys from mixing for sanitization and open tank to drain.
At EOB	Document Mix tank heel.
FOB	Once Day shift Prep has heated water loops, document Start and End times on NCCL039 for Fill machines and DW drops/hoses.

Night Lead/Tank Room Handler
Batch Change with Sanitization Standard Work

(Fill Side-All Night Shift Supervisors/Asst. Supervisors/PT's/ Handlers must be garbed and on the line at 18:55)

Time	Task to be Performed
Prior to EOB	Verify release of new batch tank, select tank, perform fill point setup and print sample tags for new batch.
At FOB	Tag sample bags to be used for SP1 (tank room samples).
At FOB	Tank Room Handler and additional handler must verify the Tank Hook Up.
At FOB	Once sanitization has been completed and sanitization filters, hoke valve caps and gray hose have been removed: <ul style="list-style-type: none"> • Perform tank valve flush • Install cap at header box • Flush SP1 and pull SP1 samples
At FOB	Once new filters are installed prior to the surge and surge tank has been drained of sanitization water then: <ul style="list-style-type: none"> • Install surge tank measuring rod • Perform flush through Filter #1 / Filter #2 and document • Perform appropriate surge tank flush per Filling 3H and drain
At FOB	Verify Filters #3 and #4 are installed downstream of surge tank.
At FOB	Begin filling surge tank with solution for FOB flush.
At FOB	Notify Machine Handlers to install flush manifolds.
At FOB	Once surge is filled flush through Filter #3 / Filter #4 and document.
At FOB	Once surge tank has filled up then signal line to begin FOB flush per Fill 4 Addendum.
At FOB	Assure entire line has flushed the required amount of solution per Filling 4 Addendum.

LPF 5, 6 and 7 Handler Batch Change with Sanitization Standard Work
 (Fill Side-All Night Shift Supervisors/Asst. Supervisors/PT's/ Handlers must be garbed and on the line at 18:55)

	Time	Task to be Performed
DAY	Prior to EOB	Have manifolds at machines, inspect for missing/ damaged tubing and correct (prior to 18:00).
DAY	Prior to EOB	Get sample cart. Stage at machine (prior to 18:00).
DAY	Prior to EOB	Get new batch type box and stage at the machine (prior to 18:00) along with die gloves and charged drill with extension/chuck. Attain a die box key.
DAY	Prior to EOB	Get BOM, Barcode setup sheet, Barcode scanner and First Imprints from the batch that is ending for imprint verification.
DAY	At EOB	Clear out machine when notified. To clear out machine: <ul style="list-style-type: none"> • Disable accumulator. • Clear machine button. • Once last pattern passes fill station then re-enable the accumulator to allow another bag pattern to load. • Once this pattern reaches the fill station (all other units in the machine should have been offloaded), then cycle stop and e-stop. • Assure only 1 complete pattern is pulled at a time.
DAY	At EOB	Sign last imprints.
DAY	At EOB	Place signed imprints (Last and First), BOM and Barcode setup sheet in red batch folder.
DAY	At EOB	Remove old type and install new type.
DAY	FOB	Change over Fill/Pack to new batch.
DAY	FOB	Verify and Document operator zone clearance.
DAY	FOB	Verify previous batch SPC is closed out.
DAY	FOB	Perform Line clearance for new batch.
DAY	FOB	Install manifold for sanitization. Crack open bleed valve on fill manifold. Begin sanitization when signaled by Tank Room Handler.
Night	FOB	At the completion of sanitization, remove flush manifold.
Night	FOB	Pull First Imprints. Sign Imprint and Imprint verification form.
Night	FOB	After acceptable first imprint is obtained or when Tank Room Handler announces-install manifold on machine and crack open drain valve on manifold, verify air inject lines removed.
Night	FOB	Tank Room Handler will signal to FOB flush per Filling 4 Addendum.
Night	FOB	After flush complete close hoke valve, remove manifold and reconnect air injector lines.
Night	FOB	If imprints not completed at this time continue with First Imprint process.
Night	FOB	If Imprint process complete and entire line has completed FOB flushes then perform FOB sampling.

RPF 1 and 2 Handler Batch Change with Sanitization Standard Work

(Fill Side-All Night Shift Supervisors/Asst. Supervisors/PT's/ Handlers must be garbed and on the line at 18:55)

	Time	Task to be Performed
DAY	Prior to EOB	Have manifolds at machines, inspect for missing/ damaged parts and correct (prior to 18:00).
DAY	Prior to EOB	Get sample cart. Stage at machine (prior to 18:00).
DAY	Prior to EOB	Get new batch type box and stage at the machine (prior to 18:00) along with die gloves and charged drill with extension/chuck. Attain a die box key.
DAY	Prior to EOB	Get BOM, Barcode setup sheet, Barcode scanner and 'First Imprints' from the batch that is ending for imprint verification.
DAY	At EOB	Clear out machine when notified by Tank Room Handler <ul style="list-style-type: none"> Assure only 1 complete pattern is pulled at a time.
DAY	At EOB	Sign last imprints.
DAY	At EOB	Place signed imprints (Last and First), BOM and Barcode setup sheet in red batch folder.
DAY	FOB	Remove old type and install new type.
DAY	FOB	Change over Fill/Pack to new batch.
DAY	FOB	Verify and Document operator zone clearance.
DAY	FOB	Verify previous batch SPC is closed out.
DAY	FOB	Perform Line clearance for new batch.
DAY	FOB	Install manifold for sanitization. Crack open bleed valve on fill manifold. Begin sanitization when signaled by Tank Room Handler.
Night	FOB	At the completion of sanitization, remove flush manifold.
Night	FOB	Pull First Imprints. Sign Imprint and Imprint verification form.
Night	FOB	After acceptable first imprint is obtained or when Tank Room Handler announces- install manifolds on machine and crack open drain valve on fill manifold, verify air inject lines removed.
Night	FOB	Tank Room Handler will signal to FOB flush per Filling 4 Addendum.
Night	FOB	After flush complete close hoke valve, remove manifolds and reconnect air injector lines.
Night	FOB	If imprints not completed at this time continue with First Imprint process.
Night	FOB	If Imprint process complete and entire line has completed FOB flushes then perform FOB sampling.

DAY Shift LPF/RPF Operator**Batch Change with Sanitization Standard Work**

Time	Task to be Performed
At EOB	Perform zone clearance of assigned machine/area. Remove any units and place in discard bin and remove batch board sheets and discard into trash or place in batch folder.
At EOB	Cut all units in the discard bin. Communicate the number of units cut to the material handler for documentation.
At EOB	Remove cut scrap/trash and place out the door to pouch side at A or E belt(Once door bar in place, door cannot be opened to pass out scrap).
At EOB	Verify handler has completed imprints then remove both the barcode foil and the black foil from the machine/printer.
At EOB	Remove air inject lines.

DAY Shift Manual Belt Operators**Batch Change with Sanitization Standard Work**

Prior to EOB	Once printed bag count is achieved clear out printer. Leave last pattern of printed bags for handler's imprints. Advance and remove the black and white used foil. Perform clearance of printer area removed all printed bags and used foil from the area.
At EOB	Remove cyclohex wick applicators from nozzle area and place on rack. (Leave 1 applicator for sample nozzle).
At EOB	Disconnect air injector lines.
At EOB	Perform clearance of nozzle area and remove all trash and bag scrap through E belt door.

Night Manual Belt Handler **Batch Change with Sanitization Standard Work**

(Fill Side-All Night Shift Supervisors/Asst. Supervisors/PT's/ Handlers must be garbed and on the line at 18:55)

Time	Task to be Performed
Prior to EOB	Get all material partial and cyclohex information and document into fill/pack materials.
At FOB	Perform FOB flush per Filling 4 Addendum when signaled by Tank Room Handler.
At FOB	Printing on Multi-up may begin-Notify operators to print.
At FOB	Pull FOB samples from manual nozzles.

Manual Belt Handler Batch Change with Sanitization Standard Work

(Fill Side-All Night Shift Supervisors/Asst. Supervisors/PT's/ Handlers must be garbed and on the line at 18:55)

	Time	Task to be Performed
DAY	Prior to EOB	Determine Bag lead needed.
DAY	Prior to EOB	Get sample cart. Stage at machine (prior to 1800).
DAY	Prior to EOB	Get new batch type box and stage at the machine (prior to 1800) along with die gloves and charged drill with extension/chuck. Attain a die box key.
DAY	Prior to EOB	Get BOM, Barcode setup sheet, Barcode scanner and 'First Imprints' from the batch that is ending for imprint verification (prior to 1800).
DAY	Prior to EOB	Sign last imprints.
DAY	Prior to EOB	Place signed imprints (Last and First), BOM and Barcode setup sheet in red batch folder.
DAY	Prior to EOB	Turn on Batch Change Light.
DAY	Prior to EOB	Chain off printer area.
DAY	Prior to EOB	Remove old type and install new type.
DAY	Prior to EOB	Clear 2199 printer area.
DAY	Prior to EOB	Document clearance and Print area operator zone clearance.
DAY	Prior to EOB	Pull first imprints from 2199. Verify/ sign imprint and imprint verification sheet.
DAY	Prior to EOB	Once all production has ended pull EOB samples.
DAY	At EOB	Verify imprints present and all documents in batch folder. Place batch fold on F belt for pouching side to pick up.
DAY	FOB	Change over Fill/Pack to new batch.
DAY	FOB	Verify and Document operator zone clearance.
DAY	FOB	Verify previous batch SPC is closed out.
DAY	FOB	Perform Line clearance for new batch.
Night	Prior to EOB	Get all material partial and cyclohex information and document materials into new batch fill/pack.
Night	At FOB	Perform FOB flush per Filling 4 Addendum when signaled by Tank Room Handler.
Night	At FOB	Printing on Multi-up may begin.
Night	At FOB	Pull FOB samples from manual nozzles.

Tank Room Prep Batch Change with Sanitization Standard Work


(Prep-All Night Shift Supervisors/Asst. Supervisors/PT's/ Handlers must be garbed and on the line at 18:50. Begins when last tank is unhooked and running out of surge)

	Time	Task to be Performed
DAY	Prior to Tank Unhook	Have the following items staged and ready following tank unhook: <ul style="list-style-type: none"> • Have sanitization filters DW flushed with whirl packs on ends. • Have 3 caps with hoke valves. • Connect gray hose to DW drop. • Disconnect filters for new batch from DW and place whirl packs on the inlets and outlets. • Obtain NCFL029 form from the fill line to document ambient flush.
DAY	After Tank Unhook	After tank is unhooked disconnect Mix Tank 21 and connect to gray hose.
DAY	Prior to Sanitization	Disconnect pipe from Mix Tank 22 and place cap with hoke valve on end directed at drain.
DAY	After Tank Unhook	Remove Filter #1 and Filter #2. Place whirl pack on the inlet piping to the surge tank. Install sanitization pipe on inlet side only.
DAY	Prior to Sanitization	Install caps with hoke valves on header box openings.
DAY	Prior to Sanitization	Tank Room Handler will shut off pumps, measure surge tank and drain. Remaining filters may be removed at this time.
DAY	Prior to Sanitization	Complete the connection between the supply pump and the surge tank.
DAY	Prior to Sanitization	Begin 6 min. ambient flush into the surge tank with the drain valve open (Document Start time on NCFL029). Bump pump on/off 2 times. At the completion document end time.
DAY	Prior to Sanitization	While performing the ambient flush install sanitization filters #3 and #4 on outlet side only with inlet side covered by whirl packs. At the completion of the ambient flush begin heating the loop and complete connections.
DAY	During Sanitization	After all connections completed begin filling surge to overflow with hot water. 30 min sanitization with filters in line required (Tank Room Handler will document).
DAY	After Sanitization	Removed filters will be tested by Night shift after water has cooled.
DAY	During Sanitization	Complete 30 minute DW sanitization of tank room drops for midweek requirement.
DAY	After Sanitization	Assure sanitization has been completed by pouch side, fill side and tank room drops prior turning loop back to cool.
Night	After Sanitization	After sanitization open surge tank to drain and remove the following: <ul style="list-style-type: none"> • Hoke valve caps • Gray hose • Piping between the supply pump and the surge tank • Filters used for sanitization
Night	After Tank Unhook	Reconnect pipes to Mix tanks.
Night	At EOB	Install new Filter #1 and Filter #2.
Night	At FOB	Install Filter #3 and Filter #4 on outlet side (inlet still covered with whirl packs), then complete the connection for the inlet side once the surge tank flush has been completed and drained by the Night shift Tank Room Handler.
Night	At FOB	Test Filters from sanitization and the previous batch.

Lines 10 / 11 Filling Tool Crib




LINE 10 – DRAWER #1



For Information Only-May Be Discarded
1 Copy Posted on Line 10


Document Owner: Nikki Buchanan

Line 10 – Drawer #2




For Information Only-May Be Discarded
1 Copy Posted on Line 10


Document Owner: Nikki Buchanan



For Information Only-May Be Discarded Document Owner: Nikki Buchanan
1 Copy Posted on Line 10




For Information Only-May Be Discarded Document Owner: Nikki Buchanan
1 Copy Posted on Line 10



For Information Only-May Be Discarded
1 Copy Posted on Line 10

Document Owner: Nikki Buchanan



For Information Only-May Be Discarded


Document Owner: Nikki Buchanan

1 Copy Posted on Line 10

Ref: L10/11 Spec. Filling Handler

Module 3

Effective Date: 08/18/17




For Information Only-May Be Discarded
1 Copy Posted on Line 10

Document Owner: Nikki Buchanan

Ref: L10/11 Spec. Filling Handler


Module 3

Effective Date: 08/18/17




For Information Only-May Be Discarded
1 Copy Posted on Line 11

Document Owner: Nikki Buchanan




For Information Only-May Be Discarded
1 Copy Posted on Line 11

Document Owner: Nikki Buchanan




For Information Only-May Be Discarded
1 Copy Posted on Line 11

Document Owner: Nikki Buchanan




For Information Only-May Be Discarded
1 Copy Posted on Line 11

Document Owner: Nikki Buchanan



For Information Only-May Be Discarded
1 Copy Posted on Line 11

Document Owner: Nikki Buchanan



For Information Only-May Be Discarded

Document Owner: Nikki Buchanan

1 Copy Posted on Line 11



For Information Only-May Be Discarded

Document Owner: Nikki Buchanan

1 Copy Posted on Line 11

Tool Re order Sheet -• **Drawer 1**

65300077- Wire Stripper

65300081- Snips

65300082- 7/16" Socket

65300083- 5/8" Socket

65300084- 1/2" Socket

65300128- 9/16" Socket

65300129- 3/4" Socket

65300130- 11/16" Socket

65300132- Ratchet

65300134- Extension

78901053- Tube Cutter

86893209- Spare Blade

86893573- V-Bar

• **Drawer 2**

65300048- 1/16" SS ALLEN WRENCH

65300049- 5/64" SS ALLEN WRENCH

65300050- 3/32" SS ALLEN WRENCH

65300051- 7/64" SS ALLEN WRENCH

65300052- 1/8" SS ALLEN WRENCH

65300053- 9/64" SS ALLEN WRENCH

65300054- 5/32" SS ALLEN WRENCH

65300055- 3/16" SS ALLEN WRENCH

65300056- 7/32" SS ALLEN WRENCH

65300057- 1/4" SS ALLEN WRENCH

65300058- 9/32" SS ALLEN WRENCH

• **Drawer 3**

65300122- 8MM Wrench

65300123- 9MM Wrench

65300124- 10MM Wrench

65300125- 11MM Wrench

65300126- 12MM Wrench

65300127- 13MM Wrench

65300106- 14MM Wrench

65300107- 15MM Wrench

65300108- 16MM Wrench

65300109- 17MM Wrench

65300110- 18MM Wrench

86893571- Micron Screen Check

• **Drawer 4**

78900101- 180° Snap Ring

78900102- 90° Snap Ring

65300085- Small Phillips

65300086- Large Phillips

65300087- Small Flat-head

65300088- Large Flat-head

65300076- Large Channel Locks

65300133- Small Channel Locks

• **Drawer 5 & 6**

65300089- 1" Wrench

65300090- 1-1/4" Wrench

65300091- 1-1/8" Wrench

Quantity:**Quantity:****Quantity:****Approval: Line 10/11 Supv:** _____

Tool Shift Verification

Week: _____

	Monday	Tuesday	Wednesday	Thursday	Friday	Saturday	Sunday
Day							
Night							

****Sign in Each block at end of shift to verify all tools are in proper location**

For Information Only

Return Completed Form to Tim Marini – May Be Discarded After Review

Line 10 Parts Kanban Standard Work

The parts Kanban has been provided to allow easier access to replacement parts that are the most commonly used on LPFs. This will enable Mechanics and PT/ Handlers to avoid downtime associated with searching for parts. This Kanban has been clearly labeled with part names and numbers and organized such that parts can be found, used, and replaced easily. Handler has been assigned days of the week that they are responsible for checking the Kanban parts cabinet and restocking. It is the responsibility of that PT/SR handler to check for needed parts and replace as needed.

- **Saturdays:** 5th and 7th will check and replace parts as needed.
- **Wednesdays:** 6th and 8th will check and replace parts as needed.
- All Shifts are responsible to receive parts from South Maintenance

For Information Only-May Be Discarded

PARENT: FILL LINES 10 & 11 SPECIFIC

PAGE 37 OF 49

Document Owner: Cathy Holloway

FL10.11 SPECIFIC8

EFFECTIVE DATE: 04/24/18

Print Fcyg:

1 Copy Posted in L10 Material Exchange Room

Go/No Go Gauge for Manual Machines Enhanced Closure Tubing

On Fill Lines 10 and 11, there are LPFs, RPFs, and Manual Fill Machines. On these machines, there are enhanced closure feeder bowls and tubing runs from the feeder bowl to the plugging station. The manual machines use tubing stock #655-00-107 which has an inner diameter of 0.380". The LPFs and RPFs use tubing stock #655-00-080 with an inner diameter of 0.375". Because the ID for both types of tubing is so similar and both tubes are clear in color, it is very easy to mix up the two. A go/no go gauge has been made to improve this problem.




The gauge is labeled with the "Manual EC Tubing"

How to Use:


Whenever you are replacing tubing on an enhanced closure feeder bowl, select what you believe to be the correct tubing from the material exchange room.

Place the gauge in to the enhanced closure tubing.

- If the gauge fits easily, you have the manual machine enhanced closure tubing.



- If you have to struggle to fit the gauge in to the tubing, or it will not fit at all, you have LPF/RPF enhanced closure tubing.



If your first tubing choice doesn't match the machine that you want to put it on, test another roll of tubing until you find the correct one.

*Note: If this gauge is ever lost or damaged, it will be set up as a stock item and can be found in the stock room.

Lockout for Solution Pump Standard Work

Description: To allow the Solution Transmission System to be safely opened when a break in fill line occurs behind the shut off valve or the break affects the shut off valve itself.

This is not a lockout procedure to be used for making any kind of electrical repairs. Lockout /Tagout per plant safety procedures should be followed and lockout/ tagout should be performed by authorized individuals only.

When to Lockout Pump: The pump should be locked out anytime solution cannot be shut off with a divert valve or cutoff valve that is a part of the Solution Transmission System or for pump repairs.

Why: Solution pump lockout is performed to prevent the solution pump from being accidentally turned back on while the solution transmission system is open.

How:

1. Stop the solution pump by pressing the red stop button located on the control panel that houses the tank controls.
2. Place the switch located at the pump plug receptacle to the **OFF** position.
3. Once the switch is in the OFF position, place a lock on the switch assuring the lock is closed completely. (Locks can be attained from the Maintenance Department).
4. The solution line can now be safely opened. Always follow correct procedures in the applicable Standard Operating Procedure for the type of break in fill line that is performed. Reference the applicable SOP for solution flushes and sampling required.
5. Once repairs are completed and the solution line is closed, the person who placed the lock on the switch may now remove the lock and turn the switch back to the **ON** position and restart the solution pump.

Standard Work for Rework % Reduction

1. Record at start of shift the filling, plugging, and rework percentages on LPF/RPF Report Worksheet
2. Continue to record data approximately every 2 hours (i.e. Quality Check).
3. If any filling, plugging, and/or rework percentages are not moving toward the goal of filling / plugging at or above 99% and rework at 3.00%, notify Line Supervisor and/or Maintenance immediately.
4. If percentages continue to not meet goal notify Line Supervisor and they will notify Maintenance Supervisor and Manufacturing Superintendent.
5. During the shift, if the filling, plugging, and/or rework percentage does not drive a favorable movement toward goal the issue will be escalated to the Shift Superintendent/Focus Force Team.

Cyclohex Dispensing Standard Work



PPE
Available
At Point Of
Use



Employee
wearing PPE and
properly
dispensing
cyclohex



Pump
Drum
Before
Dispensing

FOR INFORMATION ONLY-MAY BE DISCARDED
1 COPY POSTED IN SOUTH GARB ROOM

DOCUMENT OWNER: TIM MARINI

Line 11 Piping Verification Talley Sheet

Mix Tank to Surge Tank	A-B	B-C			
Totals					
Surge Tank to Filter Leg 1 (P1)	I-Z (3')	Z-(P1 In) (3')	(P1 Out)-N (2')	R-S (3')	S-X (3')
Totals					
Surge Tank to Filter Leg 2 (P1)	I-Z (3')	Z-(P1 In) (3')	(P1 Out)-O (2')	P-W (3')	
Totals					
Filter Leg 1: LPF5, LPF7 & Manual Filler E-belt	LPF5 U-LPF5 (2')				

Print Fig:					
Total					
	LPF7 S-Y (2')	Y- LPF7 (1.5')			
Total					
	Manual Belt X-Filler E (1.5')				
Total					
Filter Leg 2: LPF6, RPF2 & RPF1	LPF6 Q-V (2')	V-LPF6 (1.5')			
Total					
	RPF2 T-RPF2 (2')				
Total					
	RPF1 W-RPF1 (2')				
Total					

LINE 10 FILLSIDE – 6S DAILY SCORECARD

Fill Date _____ Shift _____

	YES	NO
1. All alcohol bottles are labeled and contain red color code		
2. All scrap containers are labeled and in designated area per diagram.		
3. All empty totes and lids are stored in designated area.		
4. Water hose is in designated area, hung and drained properly.		
5. All quality check gauges and scanners are in designated area.		
6. Floor is clear of dropped bags, barcode/black foil, paper towels, etc.		
7. All LPF squeegees are present and in designated area.		
8. All desks are neat and only contain line specific items.		
9. All buggies and carts are in designated area.		
10. Bag room is orderly and all items are stored in designated areas.		
11. Check all doors for free operation. Check all door pins (hinges) for signs of wear, rust and debris.		
Score: Goal 90% (# of Yes Items ÷ 10 X 100)		

Document Owner: T. MARINI

For Information Only – May Be Discarded

3 Copies Posted on Fill Side Line 10

LINE 11 FILLSIDE – 6S DAILY SCORECARD

To achieve 6S Level “2” we must comply with each of the following items.	YES	NO
1. All scrap containers are labeled and in their designated area.		
2. All water hoses are labeled, disconnected and hung on hose rack. (N/A at Batch Change)		
3. All desks, Machine HMI's and LPF cleaning stations contain only designated items and are free of clutter.		
4. All safety guards are in place.		
5. All catch pans/bins are in place and drain hoses are directed to drain.		
6. All alcohol bottles are labeled and contain the red color.		
7. All buggies, carts, ladders are in their designated areas.		
8. Line tool box and setup boxes contain only designated items.		
9. Floor, tops of equipment and fillstation cleaning tray (LPFs) is clear of trash, excessive closures and broken parts.		

Score = Pass/ Fail

Document Owner: N. Buchanan For Information Only – May Be Discard

FOB/Shift checks

All quality checks have been performed and entered for period:

1) 0700-0900

2) 0900-1100

3) 1100-0100

4) 0100-0300

5) 0300-0500

6) 0500-0700

FOR INFORMATION ONLY – MAY BE DISCARDED AFTER USE
DOCUMENT OWNER: SHEILA JONES

LPF and RPF

Maintenance, Flush and Sampling Notification

During daily production unscheduled maintenance activities or scheduled preventative maintenance may occur. It is the responsibility of the production handler/ production tech. (PT) to assure that the appropriate machine flush, sampling and follow-up testing is performed per appropriate procedure in each instance. If maintenance is performed by a mechanic then the production handler or PT should be notified of the activities prior to restart of the equipment. To aide in this process Communication Placards must be placed at the load stations of each affected piece of equipment when these activities occur. The following will give a general outline for placard use. Any additional situations that may occur should directed to the line supervisor or delegate.

- **40 Micron Mesh Screen Preventive Maintenance (PM) is scheduled to be performed:** Once production has finished on the piece of equipment the placard “Stop. BFL Flush and Samples Required” should be placed at the designated location at the load station.
- **Monthly Machine Maintnenace (PM) is scheduled to be performed:** Once production has finished on the piece of equipment the placard “Stop. Maintenance of Equipment Performed. Follow-up Testing Required”.
- **Unscheduled maintenance (emergency maintenance) or equipment adjustments:** At the time is equipment is cycle stopped or cleared of units the placard “Stop. Maintenance of Equipment Performed. Follow-up Testing Required”.

At the completion of machine flush, sampling and/or follow-up testing the placard will be removed and placed at the designated placard storage area for each machine.